ES

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### APPLICATION OPTIONS

The heating pad is designed for use on the following body parts: arms, legs, stomach, shoulder, neck, head, etc.

### GENERAL PRECAUTIONS

- After removing the packaging, check the appliance to make sure that there are no signs of damage caused during transportation. If in doubt, do not use the appliance and take it to an authorised Service Centre.

#### Risk of electric shock!

- Connect and operate the appliance only in accordance with the specifications on the rating plate.
- Examine the appliance frequently for signs of wear or damage. If there are such signs, if the appliance has been misused or does not work, return it to the supplier before switching it on again.

#### Do not use adapters, multiple sockets and/or extension leads.

- If the plug on the appliance and the mains socket are incompatible, get a qualified person to replace the plug/fuse suitable one.

#### 5. The heating pad is not to be used by children or sick children who are insensitive to heat. The electric heating pad should not be used unless the control has been switched on and set by the parent or guardian or the child has been shown how to use the device correctly.

- Do not use disabled persons.

#### To prevent injury, the appliance may be repaired by our customer service only.

#### DO NOT:

- bring the power cord into contact with hot parts.
- pull the power cord over sharp edges.
- allow heating pads in dry rooms only (not in the bathroom, etc.).
- Never use a damp or wet heating pad; wait until it has dried completely.

#### 7. This appliance must be used exclusively as an electric heating pad, as stated in the instructions. Any other use will be considered hazardous. The manufacturer is not responsible for any damage resulting from inappropriate or incorrect use.

#### IMPORTANT: Under no circumstances can the inside part of the appliance be washed or dry cleaned.

### SWITCHING ON THE APPLIANCE

The crossed wheel bin symbol shown on the appliance indicates that when it comes to dispose of the product it must not be included in household refuse. It must be taken to a special collection point for electrical and electronic appliances or returned to the distributor when purchasing a similar appliance. Under current refuse disposal regulations, users who fail to take discarded appliances to special refuse collection points may be penalised. Contact your local refuse service or the shop where the product was purchased. Your local town council can provide you with information about how to dispose of obsolete appliances.

#### B1

- When the temperature is set on the thermostat (3), the heating pad switches on and the indicator light (4) comes on.

#### Cautions:

- If the thermostat (3) is set to position "0", the heating pad switches off.
- When the heating pad ensures that a little air can circulate between the heating pad and skin. Never use the heating pad without its fabric cover.

#### Automatic switch-off

- To prevent unsupervised operation, the heating pad automatically switches off after approx. 90 minutes continuous operation. The indicator light comes on again. To switch on the heating pad again, first set the thermostat to "0" and then to the required temperature.

#### 6. Do not use pull the mains plug.

#### CLEANING, STORAGE

#### DO NOT USE THE ELECTRIC HEATING PAD IF IT IS STILL DAMP

- The thermostat must not become wet. Never hold under water. Moisture damages the thermostat.
- Do not use detergents or solvents. These may damage the appliance.

#### 8. When it comes to use any type of electrical appliance there are some basic rules that should be observed, namely:

- Do not touch the appliance with wet or damp hands or feet.
- Do not tug on the mains cable or the appliance to pull it out from the mains socket.
- Do not expose the appliance to weather conditions (rain, sun, etc.).

#### 9. Always unplug the appliance from the mains supply before carrying out any cleaning or maintenance operation on it.

#### 10. To avoid dangerous overheating, it is advisable to unwind the mains cable fully and unplug the appliance from the mains socket when not in use.

#### 11. Do not attempt to change the mains cable on this appliance. If it is damaged, the appliance must be sent away for repair.

#### 12. Hold the handle of the garment firmly when removing it from the skin. Do not pull the garment along the mains cable or use the cable as carry handle.

#### 13. If the appliance is going to be used continuously for long periods the control must be set to "0".

#### 14. ONLY to be used on the back. DO NOT lie on the heating pad as this causes it to overheat.

#### 15. Never immerse this appliance in water.

#### 16. Should you decide to dispose of the appliance it must be rendered unusable by cutting the mains cable off, after unplugging it from the mains socket. We advise you to do the same with all potentially hazardous parts on the appliance, especially those that children might be tempted to play with.

### ADVICE ON DISPOSAL:

Our goods come in optimised packaging. This basically consists in using non-contaminating materials which should be handed over to the local waste disposal service as secondary raw material.

This product complies with EU Directive 2012/19/EU.

### PICTOGRAMAS

- [No wash]
- [Wash at 40 °C]
- [Do not dry]
- [Do not iron]
- [Iron at low temperature]
- [Bleaching not permitted]
- [Do not dry-clean]

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