SIEMENS

Long Life recirculation kit LZ10AFR00



- ✓ For recirculation operation
 ✓ Including multiple regenerative active carbon filters (2 units) and a diffusor for chimless operation
- \checkmark Can be used multiple times, the regenerative active carbon filter is regenerated in the oven with a recommendation of every 4 months. Simply place in the oven for 2 hours at 200°C.

Technical Data

Dimensions of the packed product (HxWxD): ... 270 x 350 x 450 mm Standard number of units per pallet:24 EAN code:4242003774366



SIEMENS

Long Life recirculation kit LZ10AFR00

Optional accessories • For recirculation operation • Including multiple regenerative active carbon filters (2 units) and a diffusor for chimless operation Can be used multiple times, the regenerative active carbon filter is regenerated in the oven with a recommendation of every 4 months. Simply place in the oven for 2 hours at 200°C.